CURASEPT 1%
FIRST BECAUSE WE CARE

EN - Description of function

CURASEPT® Periodontal Gel is a topical gum treatment that works on a specific area and does not cause pain.

It contains PVP/VA copolymer which forms a protective film on gums that helps – by a mechanical action – to alleviate the symptoms caused by irritation to the mouth's mucous membrane.

By means of the same mechanism, it also prevents the formation and depositing of dental plaque. It includes 1% chlorhexidine digluconate which also ensures an effective antiseptic action on the film forming component.
Since there is no alcohol contained, which could cause irritant or painful effects, and with the inclusion of the A.D.S.® System, this product has a high level of acceptability.

A.D.S. ANTI DISCOLORATION SYSTEM®

A consequence of prolonged use of products containing Chlorhexidine can be discoloration of the surface of teeth.
The A.D.S.® System contained in CURASEPT® products impedes two of the main mechanisms causing this staining and helps to reduce its occurrence and extent.

Function

Topical treatment for gums to relieve soreness during dental therapy, for wearers of dentures or orthodontic braces, and for persons with periodontal disease or with implants.
It also helps prevent the formation and depositing of dental plaque.

How to use

Apply a small quantity of gel, using your finger or Curasept® Soft toothbrush, onto the gums and massage in.

To apply between teeth, use a Curasept® Specialist Mono Tuft inter-dental brush covered in gel.

In persons with periodontal disease, it is recommended to apply the product in the pockets.
If sensitisation occurs, seek advice from your dentist.

Warning

Keep out of children's reach.
Do not swallow.
Do not use this product after the expiry date shown on the tube or if the tube is not intact.
Keep at room temperature.
Please dispose of the box and its contents carefully.
Wash hands before and after use.

Contents of the box

30 ml tube ℮

CE 0373 MD

22 PAP

CURASEPT S.p.A.
Via G. Parini, 19
21047 Saronno (VA) Italy
Tel. +39 02 96 227 99
Fax +39 02 96 709 243
info@curaseptspa.it
www.curasept.it

XDP0009